TANIM

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	Manufacturer		Reference Number		Caution		Keep away from sunlight
	Date of Production		Consult instructions for use		Use by date		Keep dry
	Lot Number		Medical device		CE Certificate Notified Body Number		The product can be sterilized by steam autoclave at 134°C.
	Number of products in the box		Sterilized using irradiation		For dental professional use only.		
	Asorti Blister Package		Made of Nickel Titanium		Do not use if package is damaged and consult instructions for use		

DIRECTIONS FOR USE

PROENDO INSTRUMENTS FOR ENDODONTIC TREATMENT:

PROENDO Shaping Files (PSX, PS1, PS2)

PROENDO Finishing Files (PF1, PF2, PF3)

COMPOSITION

The cutting part of these instruments is made of a nickel-titanium alloy.

1. INDICATIONS FOR USE

These instruments are to be used only in a clinical or hospital environment, by qualified users. Application field: for the removal of dentin and shaping of the root canal.

2. CONTRAINDICATIONS

As with all mechanically driven root canal instruments, PROENDO files should not be used in cases of severe and sudden apical curvatures due to heightened risk of separation.

3. WARNINGS

This product contains nickel and should not be used for individuals with known allergic sensitivity to this material.

4. PRECAUTIONS

- Straight-line access is a prerequisite for proper root canal treatment, PROENDO files are no exception.
- Multiple use disinfection and re-sterilization cycles may lead to increased risk of file separation.
- These instruments should not be immersed in a sodium hypochlorite solution.
- Instrument reprocessing: follow the reprocessing instructions on part 7. Clean the flutes frequently during instrumentation, inspecting for signs of distortion or wear, such as uneven flutes, dull spots.
- Frequently irrigate, recapitulate and irrigate the canal throughout the procedure, minimally after using each file.
- PROENDO files should only be used in regions of the canal that have a confirmed and reproducible glide path. Establish a reproducible glide path using hand files, at least an ISO 015 size.
- Use the Shaping Files (PS1, PS2 and PSX) with a brushing action on the withdrawal stroke in order to create straight line radicular access.
- Use the Finishing Files (PF1, PF2 and PF3) with no brushing action.
- Use the appropriate finishing files to passively follow the canal to the working length, and then withdraw immediately.
- PROENDO files are manufactured with a process that results in a file that has a Gold appearance. Due to this proprietary processing, PROENDO files may appear slightly curved. This is not a manufacturing defect. While the file can be easily straightened using only your fingers, it is not necessary to straighten the file prior to use. Once inside the canal, the PROENDO file will follow the anatomy.
- Always use minimal apical pressure. Never force the files down the canal.
- For optimal usage, torque control devices are recommended.

work these instruments until a smooth, reproducible glide path is confirmed. Alternatively, mechanized glide path files (such as ProGlider® or PathFiles®) may be used after a 10 hand file.

- In the presence of NaOCl, "float" the PS1 in the canal and passively "follow" the glide path. Before light resistance is encountered, laterally "brush" and cut dentin on the outstroke to improve straight-line access and apical progression. Always brush away from the furcation.
- Continue shaping with PS1 as described until the depth of the 15 hand file is reached.
- Use the PS2, exactly as described for the S1, until the depth of the 15 hand file is reached.
- In the presence of a viscous chelator or NaOCl, scout the apical 1/3 with 10 and 15 hand files and gently work them until they are loose at length.
- Establish working length, confirm patency and verify the presence of a smooth reproducible glide path in the apical 1/3.
- Use the PS1, with a brushing action, until working length is reached.
- Use the PS2, with a brushing action, until working length is reached.
- Reconfirm working length, irrigate, recapitulate and re-irrigate, especially in more curved canals.
- Use Finishing File PF1, in a "non-brushing" action, with each insertion deeper than the previous insertion until working length is reached. Do not leave the file at working length for longer than one second.
- Gauge the foramen with a 20 hand file. If the instrument is snug at length, the canal is shaped and ready to be obturated.
- If the 20 hand file is loose at length, proceed to the PF2 and, when necessary the PF3, PF4 and PF5, with the same non-brushing motion to working length, gauging after each Finishing file with 25, 30, 40 or 50 hand files respectively.
- If necessary, use the PSX with a brushing motion to move the coronal aspect of the canal away from furcal concavities and/or to create more coronal shape. PSX can also be used to optimally shape canals in shorter roots.
- The PROENDO sequence is the same regardless of the length, diameter or curvature of the canal.

8. DISINFECTION, CLEANING AND STERILIZATION

Reprocessing procedure for dental instruments.

FOREWORD

- Devices that are marked as "sterile" do not require any specific treatment before the first use. For all other devices not labelled "Sterile", sterilization prior first use is required according to section III – STEP-BY-STEP INSTRUCTIONS part 4 to 8 of this DFU.
- For those devices that are not labelled "single use", re-processing of the devices should be carried out as per this DFU. For hygiene and sanitary safety purposes, these instruments must be cleaned and sterilized before each re-use to prevent any contamination.
- Excluded devices: Unclip and Mooser Calcifiable plastic posts cannot be sterilized and must be disinfected by immersion NaOCl (2,5 % at least) during 5 min. at ambient temperature.

GENERAL RECOMMENDATION

The PROENDO rotary files can be used at motor speeds between 250 rpm and 350 rpm. Recommended motor settings:

File Size	Speed [rpm]	Torque [N•cm]
PROENDO PS1 & PSX	300	5,10
PROENDO PS2 & PF1	300	1,50
PROENDO PF2, PF3	300	3,10

The speed and torque settings indicated in the above table are for example only and may vary according to each user preferences and motor capabilities.

5. ADVERSE REACTIONS

As with all mechanically driven root canal instruments, PROENDO files should not be used in cases of severe and sudden apical curvatures due to heightened risk of separation.

6. STEP BY STEP INSTRUCTIONS FOR PROENDO FILES

a) Radiographic Evaluation

Review different horizontally angulated radiographs to diagnostically determine the width, length, and curvature of any given root and canal.

b) Access Preparation

Create straight-line access to the canal orifice(s) with emphasis on flaring, flattening, and finishing the internal axial walls.

7. PROENDO SHAPING TECHNIQUE

- The crown down technique is the technique of choice for rotary instruments.
- Create straight-line access to canal orifice.
- In the presence of a viscous chelator (such as Glyde® File Prep root canal conditioner) passively scout the coronal 2/3 with 10 and 15 hand files. Gently

- 1-Use only a detergent solution, with disinfecting effect, which is approved for its efficacy (VAH/DGHM-listing, CE marking, FDA approval) and in accordance with the DFU of the detergent solution manufacturer. For all metal devices, it is recommended to use anticorrosion disinfecting and cleaning agents.
 - 2-For your own safety, please wear personal protective equipment (gloves, glasses, mask).
 - 3-The user is responsible for the sterilization or disinfection of the product for the first cycle and each further usage as well as for the usage of damaged or dirty devices where applicable after sterilization.
 - 4-It is safest for the practitioner to use our devices only once. Should our devices be reused, we recommend that they should not be used more than 5 times. After each processing they should be carefully inspected before use: the appearance of defects such as deformations (bent, unwound), breakage, corrosion, loss of colour coding or marking, indicate that the devices are not able to fulfil the intended use with the required safety level and must therefore be discarded.
- For our root canal shaping instruments we recommend not to exceed the following maximum number of uses:

	Manufacturer		Reference Number		Caution		Keep away from sunlight
	Date of Production		Consult instructions for use		Use by date		Keep dry
	Lot Number		Medical device		CE Certificate Notified Body Number		The product can be sterilized by steam autoclave at 134°C.
	Number of products in the box		Sterilized using irradiation		For dental professional use only.		
	Asorti Blister Package		Made of Nickel Titanium		Do not use if package is damaged and consult instructions for use		